

product care guide



**WIELAND**



## thank you for your purchase!

This care guide contains specific information to keep your Wieland products performing well. With the proper care, combined with Wieland's durable and environmentally responsible designs, you can expect years of use out of your furniture.

If you need assistance or have any concerns please contact [sales@wielandhealthcare.com](mailto:sales@wielandhealthcare.com). For warranty information visit <https://wielandhealthcare.com/about-us/warranty/>.

### material care information

---

<a href="#">2K poly wood finish</a> .....	3
<a href="#">urethane</a> .....	3
<a href="#">thermoform</a> .....	4
<a href="#">resin &amp; plastic</a> .....	4
<a href="#">solid surface</a> .....	4
<a href="#">powder coated metal</a> .....	5
<a href="#">laminat</a> e.....	5
<a href="#">fabric</a> .....	6



## 2k poly wood finish

---

Wieland's 2K Polyurethane wood finishing system is applied with the latest in computer controlled equipment. This 2-component, or 2K, polyurethane sealer and top coat has all of the advantages of traditional catalyzed varnish with the added benefits of increased durability and performance.



To test this, Wieland created an environment where each of the following noted chemicals were applied and rubbed dry 36,400 times. Each of the following chemicals were diluted to the manufacturers recommendation: Bleach (e.g. Clorox®), Isopropyl Alcohol (e.g. CleanPro®), Hydrogen Peroxide (e.g. Oxycide™), Quaternary Ammonium (e.g. Virex®) and Phenolic (e.g. Wex-cide).

After this very rigorous testing, the samples did not show a noticeable difference in Gloss Level, Film Thickness or Cross Hatch Adhesion from the sample prior to testing. This is a real testament to the superior hard use resistance inherent with this finish chemistry.

### **general cleaning instructions**

Everyday cleaning may be done with mild detergent, or soap and warm water, drying surfaces immediately after with a soft cloth. Following application of Sanitizing or Germicidal cleaning products, it is recommended to wipe the wood surface with a water-dampened cloth to remove any remaining residue.

## urethane

---

### **general cleaning instructions**

Clean with mild soap and water. Diluted cleaners and disinfectants may also be used, but wiped dry immediately after application. Do not leave to dry.





## thermoform

Thermoform will retain its original beauty with reasonable handling and care. Wieland's thermoform products may be cleaned with all standard healthcare cleaners if they are diluted according to the manufacturer's recommendation: Hydrogen Peroxide (e.g. Oxycide™), Quaternary Ammonium (e.g. Virex®) and Phenolic (e.g. Wex-cide).

### general cleaning instructions

To ensure that your surfaces maintain their original appearance, clean with soap and water, or any of the healthcare cleaners listed above.



## resins & plastic

Wieland's resin products may be cleaned with all standard healthcare cleaners if they are diluted according to the manufacturer's recommendation: Bleach (e.g. Clorox®), Isopropyl Alcohol (e.g. CleanPro®), Hydrogen Peroxide (e.g. Oxycide™), Quaternary Ammonium (e.g. Virex®) and Phenolic (e.g. Wex-cide).

### general cleaning instructions

Mild detergents are often more effective at removing general dirt and soil for everyday cleaning. Scuff marks from shoes or other hard surfaces may be removed using mild abrasive cleaners.

### repairs

The resin frame has a uniform color throughout. Any small nick or scratch can be sanded smooth.



## solid surface

Corian® was created for a lifetime of easy care. Wieland's solid surface products may be cleaned with all standard healthcare cleaners if they are diluted according to the manufacturer's recommendation: Hydrogen Peroxide (e.g. Oxycide™), Quaternary Ammonium (e.g. Virex®) and Phenolic (e.g. Wex-cide).

Visit <http://www.corian.com/-use-care-> for more information or call 1-800-4Corian (1-800-430-6072).





## powder coated metal

### general cleaning instructions

Light soiling may be removed with a damp cloth or household glass cleaners. Heel marks, heavy soiling and surface marring may be removed with automotive-grade polishing compounds or car wax. Do not apply abrasive cleaners, polish or lubricants to table bases or moving parts.

### repairs

Damaged paint may be repaired with touch-up paint which is available from the factory. Touch-up paint is supplied in powder form to provide an exact match with the existing color. To use the powder, add a small quantity of commercially available lacquer thinner to a small quantity of the powder, stirring with a metal stir until the desired consistency is achieved. The paint solution may then be applied using a brush. Avoid using too thick of a consistency that may result in brush marks. Allow at least one hour to dry.



## laminates

Laminate products may be cleaned with all standard healthcare cleaners if they are diluted according to the manufacturer's recommendation: Bleach (e.g. Clorox®), Isopropyl Alcohol (e.g. CleanPro®), Hydrogen Peroxide (e.g. Oxycide™), Quaternary Ammonium (e.g. Virex®) and Phenolic (e.g. Wex-cide).

### general cleaning instructions

To clean the surface, use a damp cloth or sponge and a mild soap or detergent. Rinse thoroughly with warm water and wipe dry.



## fabric

### general cleaning instructions

Please refer to the respective fabric supplier for cleaning instructions on fabrics.



**WIELAND**

888.WIELAND  
wielandhealthcare.com