



amari

WIELAND



14981
recliner

simplicity redefined

Amari's comprehensive simplicity is deliberately designed to translate style into sustainability.

Easily replaceable components and covers extend its life span in today's demanding healthcare environments. High density, contoured foam cushions offer superior, long-lasting comfort.

Cleaning is effortless with the optional removeable seat, and hard surface panel and arm cap. Amari's back reclines independently from its ottoman, resulting in multiple comfort options. Patients can select the perfect recline angle while choosing to raise or lower their feet.

Push red pedal for exceptional braking power; green to activate auto-align tracking casters; or neutral position for casters with zero-turn radius



14981
central brake

highlights

- independent back motion so orthopedic patients can recline comfortably without raising legs
- all components are easily renewed on site
- amari's small footprint efficiently reserves space for necessary equipment
- large casters provide 5" of clearance for patient lifts and mobility devices*
- moisture barrier is standard
- 500 lb. weight capacity
- 25 year warranty; 10 years on mechanism, 3 years on casters

hard surface arm panels (solid surface or thermoform) are easy to wipe down and impervious to harsh healthcare cleaners



14W981
recliner, wide


fold down arm for patient transfer in both the seated and supine position



14981
single transfer arm

*does not apply when footrest option is specified

customized comfort




Patients and visitors will easily relax as amari's back and ottoman operate independently of each other.

options and accessories



urethane arm cap
(standard)



solid surface cap



arm panel
thermoform
(standard)



arm panel
solid surface



quick release seat



full flip table
thermoform or
solid surface



half flip table
thermoform or
solid surface



foley bag hook



pushbar
urethane



footrest
black only




single transfer arm
available either side



double transfer arm



central brake



red: brake engaged
green: tracking casters activated
neutral: zero-turn casters
(as shown)

specifications

Product No.	Description	W	D	H	SW	SD	SH	AH	WT	YDG	WT CAP
14981	recliner	29	34.25	46.75	22	20.75	20	26.75	113	4	500
14W981	recliner, wide	32	34.25	46.75	25	20.75	20	26.75	125	4.5	500

888.WIELAND
wielandhealthcare.com