

# Accord™ recliner

Nightingale Awards

0

## nurse health & patient safety

Partnering with patients, clinicians and designers, the accord<sup>™</sup> recliner addresses caregiver health, patient safety and design flexibility. Wieland's innovation is apparent in multiple new-to-the-world features.

#### safeBrake

Brake pads drop down and grip the floor securely, delivering exceptional stability. An industry first, Wieland separates the casters and brakes. Accord's dedicated brake system outperforms brake systems that are combined with the casters.

- makes a solid stop on uneven floors
- greater surface area increases staying power
- brake retracts completely when not in use
- brake pads stay clean

#### intelliTrac™

Magnetic casters auto-align and auto-release offering the control to change direction quickly from straight to sideways.

- track straight or move perpendicular with a nudge
- designed to work with the safeBrake system, casters will not rotate or roll when brake is engaged
- available in black or grey

#### quick release seat

Provides easy access for cleaning and infection control.



#### transfer arm

A front-released arm that pivots towards the back of the chair for an obstacle-free, hands-free patient transfer.

#### controlCove

Centrally located foot pedals that are conveniently accessed from the front–no bending, stooping or reaching required.





sideways



straight



## highlights

- safeBrake provides exceptional stability by gripping the floor rather than merely stopping the casters
- the intelliTrac<sup>™</sup> caster system helps guide the chair straight down a corridor to minimize strain
- guick-release seat option supports infection prevention and safety efforts
- centrally located foot pedals give access to all controls from the front of the chair
- patient transfers are made easier with the shortest transfer distance between the bed and the chair
- a full length, user-activated sleep position reclineAssist bar is located neatly under the arm cap
- faceted back features a slimming two-panel design creating a waistline while giving the chair an open feel
- beveled heal guard encourages proper nose-over-toes balance during egress
- the fold-down table has been designed to stow slowly, avoiding sudden drop down movements
- the fold-down table comes standard with a spill groove preventing small items from rolling off to the floor
- high strength, 10-gauge steel recliner frame with welded arm support for rigidity and longevity
- soft edges and corners with thick, dense foam prevents damage to walls and upholstery
- 25-year warranty; recline mech, transfer arm, intelliTrac casters and safeBrake–10-year warranty; casters–3-year warranty; electrical-5-year warranty





## standard features



reclineAssist bar Sleep and Trendelenburg only



arm cap Wood, urethane, thermoform, or solid surface



oversized ottoman t-shaped seat cushion allows ottoman 30% larger surface area



arm insert panel Fabric or thermoform



faceted back design



soft seat



controlCove Black, sable, dove or silver metal powder coat





locking casters Black or grey



safeBrake



intelliTrac<sup>™</sup> with safeBrake Specify black or grey casters



transfer arm Specify right, left, or both



pushbar Black, sable or dove



retractable footrest Black only



Thermoform or solid surface;



quick release seat




headrest Specify fabric



USB-A/C Power Left-hand facing only. Available with or without Drip Guard. Specify silver or black



orthopedic height seat 21" seat height. Available with safeBrake option only



heat + massage Heat or heat + massage options available. Controls located on the inside of the right arm (as facing)



IV pole Standard





Black only



fold-down table Specify right, left or both. Specify white, dolphin, nutmeg or charcoal



## compact and comfortable

The Accord sleep chair provides a comfortable and convenient sleeping solution for overnight guests, creating a supportive atmosphere that reduces anxiety for both patients and loved ones. The Accord sleep chair fosters a healing environment that accelerates recovery, even in patient rooms with a smaller footprint.



## highlights

- contains a one-step mechanism for ease of use
- extends to a 72" sleep surface with an intuitive pull handle in one motion
- arm caps available in solid surface, urethane or thermoform
- comes standard with moisture barrier and fully renewable components
- standard moisture barrier protects foam from spills and liquids
- wide scale model available
- casters are standard
- 25-year warranty; 10-year warranty on sleep chair mechanism



specifications Product No. Description W D H SW SH AH WT YDG WT CAP SD 65981 Recliner 34.25 28.75 46.25 21.5 21 20 26.5 142 5.5 500 65W981 Recliner, Wide 34.25 46.25 23.5 21 5.75 500 30.75 20 26.5 146 65981S **Sleep Recliner** 28.75 34.25 46.25 21.5 21 20 26.5 150 5.5 500 65W981S Sleep Recliner, Wide 30.75 34.25 46.25 23.5 21 20 26.5 154 5.75 500 65981T Trendelenburg Recliner 28.75 34.25 46.25 21.5 21 20 26.5 158 5.5 500 65W981T Trendelenburg Recliner, Wide 34.25 21 30.75 46.25 23.5 20 26.5 162 5.75 500 65841 Sleep Chair 30 33.5 33.25 24 21 18 25 110 4.75 400 65W841 Sleep Chair, Wide 36 33.5 33.25 21 18 25 128 6.25 400 30





888.WIELAND wielandhealthcare.com